

# **Cross Check™**

# Torque Seal® Tamper-Proof Indicator Paste

#### **HANDLING**

- Obtain special instructions before use.
- Do not handle until all safety precautions have been readand understood.
- Do not handle, store or open near an open flame, sources of heat or sources of ignition.
- Protect material from direct sunlight. Explosion-proof

#### **DIRECTIONS**

- Knead tube before use.
- DRY TIME: Dry to touch/handle in 30-40 minutes.
- 80% curing occurs between 48-72 hours. Completion of curing takes place between 7-10 days.
- Exposing to elevated temperature OR blowing hot air can speed up the drying/curing process.

### **STORAGE**

• Store locked up. Keep away from heat, sparks and open flame. Prevent electrostatic charge build-up by using common bonding and grounding techniques. Store in a cool, dry place out of direct sunlight. Store in tightly closed container. Store in a well-ventilated place. Keep in an area equipped with sprinklers. Store away from incompatible materials (see Section 10 of the SDS).

## **DISPOSAL**

 Collect and reclaim or dispose in sealed containers at licensed waste disposal site. Dispose ofcontents/container in accordance with local/regional/national/international regulations.

### **FEATURES**

- Tamper Proof Torque Mark
- Safe for most surfaces
- Highly visible- blue, yellow, red, green, and orange colors are fluorescent under blacklight
- Excellent adhesion
- Visually detects loose nuts & bolts
- Provides warranty protection
- Meets BOEING Specification, BMS8-45 Type II



# **PHYSICAL PROPERTIES**

Appearance / Physical State	Liquid	
Odor	mild	
Flash Point	105.1 °F (40.6 °C)	
Boiling Point / Range	276.98 - 485.06 °F (136.1 - 251.7 °C	
VOC Content	42.28%, 430 g/L	
Vapor Density	> 1 (air = 1)	
Prop 65 Warning	YES	

#### **AVAILABLE IN**

Net	Cnt.	Units/	Part #	
Contents	Type	Case	US	Canada
1 oz	tube	24	83314	C83314
1 oz	tube	24	83315	
1 oz	tube	24	83316	C83316
1 oz	tube	24	83317	C83317
1 oz	tube	24	83318	C83318
1 oz	tube	24	83319	C83319
1 oz	tube	24	83320	C83320
1 oz	tube	24	83321	
	1 oz	Contents Type   1 oz tube   1 oz tube	Contents Type Case   1 oz tube 24   1 oz tube 24	Contents Type Case US   1 oz tube 24 83314   1 oz tube 24 83315   1 oz tube 24 83316   1 oz tube 24 83317   1 oz tube 24 83318   1 oz tube 24 83319   1 oz tube 24 83320

