

TECHNICAL DATA SHEET

SCREEN2P Wear® Icaridin

Insect Repellent for the Impregnation of Clothes

DESCRIPTION

Vector-borne diseases such as Yellow Fever, Zika, Malaria, Chikungunya and Leishmaniosis are killing over 1 million people every year, with just 1 bite.

With bites occurring at day as at dawn, sunset and night, and with 40% of the mosquito bites occurring through clothing, SCREEN2P Wear® Icaridin is the ideal long lasting protection needed on clothes by the travellers to stay safe during all they journey and stay. SCREEN2P Wear® Icaridin exclusive non-flammable, Icaridin based insect repellent for clothes and fabrics, provides 8 day proven efficacy.

SCREEN2P Wear® Icaridin can be used on all types of fabrics. Its pleasant water based formulation, create an efficient protective screen against the inconvenience and risks of insect bites.

FEATURES & BENEFITS

- Easy to apply exclusive Icaridine based formula.
- Alcohol free, water based, non-flammable formulation.
- Protective screen against mosquitoes (*Aedes, Anopheles and Culex*).
- Demonstrated protection for 8 days.
- Does not leave marks or stains on the treated fabrics.
- For all types of fabric: cotton / synthetic, clothing, denim or combat uniforms, curtains, nets, etc.
- Lotion that can get in contact with skin (except mucous)
- Discreet floral odour.
- Leak proof bottle thanks to crimped pump.
- Different practical format from 30mL to 100mL.
- Recyclable packaging.



For complete protection, use jointly with our SKIN2P Body® insect repellent for the skin.

APPROVALS

- Efficacy proven in accordance with protocol described by the WHO for repellents (WHO/HTM/NTD/WHOPES/2009.4), World Health Organization, 2009

CERTIFICATIONS AND ACCREDITATIONS

- PSA certified to ISO 9001 version 2015
- PSA NATO Cage Code : FA2X6

PHYSICAL PROPERTIES

Active ingredients : 20.6% w/w Icaridin n° CAS 119515-38-7
0.9% w/w Geraniol n° CAS 106-24-1

Appearance: White liquid
Odour: Slight floral odour
Shelf life : 36 months

TECHNICAL DATA SHEET

APPLICATIONS – USE

- The spray delivers 0.3mL of product (~0,3g) per pump.
- Dose to apply: 17mL of product per 1 m²
- The 100mL Screen2P Wear® Icaridin spray pump can treat 5.9 m² of fabric (equivalent to more than 4 sets of pants and shirts).
- The 30mL bottle enables the treatment of 1,8m² of fabric (equivalent to 2 shirts + 1 pants).
- For optimal and practical spraying, place the garment vertically on a hanger.
Shake well and then spray 15 cm from the garment; A single spray in the same place is enough to achieve the dose required for efficacy.
To guarantee even distribution of the product on the garment carry out juxtaposed sprays, i.e. the spraying rounds must touch each other.
Repeat the treatment on the back.
- For garments that are laid flat:
Tilt the bottle slightly and spray at a distance of 15 cm from the garment.
Impregnation is checked visually by the change in colour due to the fabric becoming humid.
- The drying time is 30mins to 1 hour, depending on the temperature.
- Renew the application after 8 days.

Test the product beforehand on a small portion of fabric. Wash hands before eating.

Use biocides with caution. Before use, read the label and information concerning the product

ORDERING INFORMATION

Product Code	Packaging	Label	Units /Carton	Carton Weight (Kg)	Carton Dimensions (cm)	CO ₂ /units (EXW Paris)
AR0000040	100mL spray	French & English	50	6,3	43 x 22 x 16	Coming soon
AR0000066	30mL spray	French & English	100	4,9	33 x 33 x 10	Coming soon

Looking for a specific format or labelling ? Please contact us.



**CONTACT US TO KNOW MORE ABOUT CARBON FOOTPRINT
AND CARBON OFFSET FOR CO₂ NEUTRAL PRODUCTS**



STORAGE

Store in a cool dry place. Do not store in direct sunlight

DISPOSAL

Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

WARRANTY – All statements, information and data presented herein are believed to be accurate and reliable but are not to be taken as a guarantee, expressed or implied, for which seller assumes legal responsibility and they are offered solely for your consideration, investigation and verification. Statements or suggestions concerning possible use of this product are made without representation or warranty that any such use is free of patent infringement and are not recommendations to infringe on any patent. Created 26 June 2019. Modified 1st March 2024. Date Printed 13/03/2024 7:46 PM